


Master the extreme



NobelPro Line™

NobelPro Line™

Master the extreme with our outstanding, most comprehensive treatment range & training. From common cases to severe resorption, you have an opportunity like no other to treat edentulous or soon to be edentulous patients. NobelPro Line addresses the most challenging cases with the All-on-4® treatment concept.

Adapt for your patients
Be flexible, be agile. Take your practice to the next level with a graftless solution featuring the most comprehensive product range for edentulous solutions. NobelPro Line allows you to respond to your patient’s unique anatomy and master even the most challenging maxillary cases.

Three steps to mastery
Building on your existing skills with the All-on-4® treatment concept, our three-stage program is carefully designed to take you to the next level of clinical competency. After completing the program we invite you to join our exclusive group, **The Circle™**, providing you with a unique opportunity to collaborate with global experts.



PRODUCTS

Expand your treatment options in moderate to severe bone resorption; choose from our extensive range of product options from NobelParallel and NobelSpeedy to NobelZygoma implants, NobelProcera implant bridges, bars and a full range of Multi-unit Abutments.



TRAINING

Expand your abilities in edentulous cases through training individually tailored for your professional needs. Through our three-stage program, gain confidence, support and skills to become prepared for your first complex cases.




MENTORING

Throughout the NobelPro Line program, a dedicated mentor will support your progression and provide guidance with real-life cases. You will also visit your mentor’s practice and learn from their expertise in action.



NobelProcera® Restorations
Fully customized prosthetic solutions



All implants designed with special apex for bicortical anchorage



Implants

NobelParallel™CC
Internal connection

P
7–18mm
Ø NP 3.75mm; RP 4.3, 5.0mm

NobelSpeedy®
External connection

S
7–25mm
Ø NP 3.3 mm; RP 4.0, 5.0mm
Short and narrow
7.0 mm
Ø NP 3.3mm
Wide Body
8.5–18.0mm
Ø RP 5.0 mm

NobelZygoma™
External connection

Z
30.0–52.5 mm
0°
30.0–50.0 mm
Ø RP 5 mm
45°
30.0–52.5 mm
Ø RP 3.9 mm

Multi-unit Abutments

Xeal™ 0° 17° 30°

0° 17° 30°

0° 45° 60°
45° 0° 17°



Take your skills in the All-on-4® treatment concept to the next level. Join the global community of NobelPro Line™ qualified professionals.

Register your interest at
nobelbiocare.com/proline

GMT 64017 © Nobel Biocare Services AG, 2019. All rights reserved.

Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.