



BLACK LINE RUBBER DAM CLAMPS

| An Exceptional Restoration Begins with Excellent Isolation

Hu-Friedy's 48 clamps, made in the USA from proprietary stainless steel, provide full-mouth coverage.

POINTS OF PERFORMANCE


- **Satin Steel®** which resists corrosion while aiding in flexibility and memory
- **A dark coating** allowing for improved intraoral contrast and visual acuity
- **Matte finish** for decreased light reflection and glare



ANTERIOR RUBBER DAM CLAMPS



MOLAR RUBBER DAM CLAMPS (GENERAL)




COMPLEMENTARY PRODUCTS



Rubber Dam
Punch-Ainsworth
RDPA

Rubber Dam
Forceps
RDF



Rubber Dam Frame, Child
RDCF5




Rubber Dam Frame, Adult
RDAF6


PREMOLAR RUBBER DAM CLAMPS



UPPER MOLAR RUBBER DAM CLAMPS



LOWER MOLAR RUBBER DAM CLAMPS



Visit us online at hu-friedy.eu

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA
Customer Care Department: Free Phone 00800 48 37 43 39 | Free Fax 00800 48 37 43 40
Website: Hu-Friedy.eu

©2018 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-5696GB/0518

