

#whdentalwerk



video.wh.com

Now at  
your dealer  
or [wh.com](http://wh.com)



# The new solution in **the Zygoma** **application**



## Unique, easy and safe

The specially for the insertion of Zygoma implants developed 20:1 handpiece facilitates the work of the surgeon. The angulation is unique and allows a much better view of the treatment site. So the surgical procedure is also possible to difficult accessible positions in the zygomatic bone.

### Easy access to treatment site

The angle of the handpiece allows optimal access for difficult interventions in the area of the zygomatic bone.



### High application security

The world's only 20:1 handpiece with contra-angle chucking system ensures application security at high concentricity.



### Fatigue-free working

The ergonomic design of the handpiece ensures long and fatigue-free working.




### Efficient cooling

The spray exits near the instrument's area of activity and cools the operating site extremely efficiently.



## Technical data Zygoma

Type	SZ-75
Transmission ratio:	20:1
Coupling system:	ISO 3964 (DIN 13940)
External cooling system:	single
Chucking system:	twist chuck system
Rotary instruments:	Zygoma burs, surgical burs and cutters with contra-angle shank Ø 2,35 mm
Maximum drive speed:	50,000 min <sup>-1</sup>
Can be dismantled:	yes
Indications, e.g.:	Implantation of Zygoma implants osteotomy





### **Dr Armando Lopes**

Specialist in Oral Surgery & Implantology  
MALO CLINIC, Portugal

“W&H offers the best handpiece for the insertion of Zygoma implants: The unique combination of the specially angled handpiece, the efficient spray cooling and a secure chucking system allows maximum comfort. For doctor and patient.”



### **Dr Ana Ferro**

Specialist in Periodontology  
MALO CLINIC, Portugal

“Beautifully machined for quality and balance. W&H sets new standards in the Zygoma application. Incredibly easy to use with Implantmed and fatigue-free working thanks to the ergonomic design of the handpiece are guaranteed.”



**W&H Dentalwerk Bürmoos GmbH**

Ignaz-Glaser-Straße 53, Postfach 1

5111 Bürmoos, **Austria**

**t** +43 6274 6236-0

**f** +43 6274 6236-55

office@wh.com

**wh.com**