


Sinus Instrument Kit




24.995.01

ZEPP Sinus Instrument Kit acc. to Dr. Stiller



Soft shank for an individual bending



Sinus Instrument Kit

acc. to Dr. Stiller

Immediate sinus floor elevation with or without bone cover

The **HELMUT ZEPP** Sinus Instrument Kit has been developed together with Dr. Stiller and is particularly suited for difficult maxillary sinus structure (septa, maxillary pillars, scarred mutations in change with intact maxillary sinus mucosa).



Curriculum Vitae PD Dr. med. Dr. med. dent. Michael Stiller

born September 4, 1961

1983-1988 Academic studies of dentistry at the Charité

1988 Residency in oral surgery, 1994 dissertation as Dr. med. dent.

1994 Start of activity as Deputy at the Clinic and Polyclinic for Dental and Oral Medicine at the Free University of Berlin, Department of Dental Surgery and Radiology

2001 Habilitation in Dental and Oral Medicine

since 2000 Oral Surgeon in private dental and medical practice, focus on implantology

since 2005 Research activities within the DFG (German Research Foundation)



HELMUT ZEPP
MEDIZINTECHNIK GMBH

NEWS

10/12



**Surgical Aspirator,
SinusLine, titanium tip,
slotted**

Distributed by:

Sinus Instrument Kit
















Sinus Instrument Kit

acc. to Dr. Stiller

Advantages of the system:

- Instruments which are adapted and adjusted perfectly to the anatomical conditions in the maxillary sinus
- Flexible working ends (instruments are pre-bended for the normal antral anatomy)
- Two different kinds of instruments with blunt and sharp edges for dissecting mucosa on plane and rough internal surface of the maxillary sinus
- Delicate working ends for most complicated situations during the operation, for example cicatrization, situation after mouth-antrum connection, antral dissection of maxillary sinus mucosa on the root surface and very thick facial maxillary sinus wall, dissection of distinctive costodiaphragmatic recess
- Working ends which do not glare on burnished edge, therefore it is possible to work atraumatic and ergonomically
- Beneficial and reduced self-swinging behavior of the instruments despite the flexibility of the working ends
- Two mixing cups to collect autologous bone and to mix blood with bone replacement material
- Special suction tubes for the maxillary sinus with perforations to avoid suction perforation; due to the delicate working ends it is possible to work accurately, also in difficult situations



	Article description	Order quantity
	24.995.01 Sinus Instrument Kit acc. to Dr. Stiller consisting of:	[]
	85.185.12 Tray with lid incl. grading for Sinus Kit, stackable	
	08.902.031HF HB-Drill, round	[]
	08.906.023C Diamond, round 023C, Ø 2.3 mm, hand piece, total length 44.5 mm	[]
	08.906.029C Diamond, round 029C, Ø 2.9 mm, hand piece, total length 44.5 mm	[]
	19.651.13 Surgical Aspirator, SinusLine, with titanium tip, Ø 1.5 mm, slotted	[]
	19.651.14 Surgical Aspirator, SinusLine, with titanium tip, Ø 3.0 mm, slotted	[]
	41.822.01 Stiller Sinus-Elevator, blunt, 4.0 mm, titanium ZEPF-Line, green, double-ended, 19.5 cm	[]
	41.822.02 Stiller Sinus-Elevator, blunt, 2.5 mm, titanium ZEPF-Line, green, double-ended, 19.5 cm	[]
	41.822.03 Stiller Sinus-Elevator, blunt, 3.0 mm, titanium ZEPF-Line, green, double-ended, 19.5 cm	[]
	41.822.04 Stiller Sinus-Elevator, universal, 3.0 mm, titanium ZEPF-Line, gold, double-ended, 19.5 cm	[]
	41.822.05 Stiller Sinus-Elevator, universal, 6.0 mm, titanium ZEPF-Line, blue, double-ended, 19.5 cm	[]
	41.822.11 Stiller Sinus-Elevator, sharp, 4.0 mm, titanium ZEPF-Line, red, double-ended, 19.5 cm	[]
	41.822.22 Stiller Sinus-Elevator, sharp, 2.5 mm, titanium ZEPF-Line, red, double-ended, 19.5 cm	[]
	41.822.33 Stiller Sinus-Elevator, sharp, 3.0 mm, titanium ZEPF-Line, red, double-ended, 19.5 cm	[]
	41.868.07 Kirsch Sinus 7 Plugger, titanium ZEPF-Line, blue, double-ended, 17.5 cm, Spoon 8x10 mm	[]
	85.251.00 Mixing Cup, stainless steel, incl. plastic lid, Ø 4 cm (2 pcs. included in the kit)	[]



HELMUT ZEPF
MEDIZINTECHNIK GMBH

Order by fax to +49 (0) 7464 / 98 88 88 or online at: www.zepf-dental.com

Customer Number

Company / Practice

Street, Number

Postal Code, City

Date

Source of Supply / Depot

Phone / Fax

Signature of the Client