

# Rescue instrumentation

## Product overview

---

### Abutment / clinical screw removal

In the event that an abutment screw breaks, there's a safe solution, even when part of the screw remains in the screw channel.

While the Abutment Screw Remover will resolve the majority of such instances efficiently, Reverse Drills and Abutment Retrieval Instruments enable retrieval in more complex cases.

---

### Implant retrieval

Managing a retrieval case where the implant driver no longer engages?


The Implant Retrieval Instrument is designed precisely for this scenario, regardless whether it's due to implant fracture, a damaged connection or strong osseointegration. For more difficult cases, Trephine Drills are available.

---

### Abutment retrieval

Removal of an abutment held by a tight conical connection is simple with the Abutment Release Pin.

For titanium and zirconia restoration, including the NobelProcera ASC Abutment and the FCZ Implant Crown, removal is simplified with a range of dedicated abutment retrieval instruments.



# Surgical components

## Implant retrieval

### Implant Retrieval Instruments

	Conical Connection	Tri-channel	External Hex
22 mm	37471	NP, RP	NP, RP
31 mm	37472	NP, RP	NP, RP
31 mm	37473 3.0		
22 mm	37474 NP		WP
22 mm	37475 RP	WP	
22 mm	37476	6.0	
22 mm	37927 WP		

STERILE R



### Implant Rescue Collar

	Tri-channel
Ø 3.5 mm	37478 NP
Ø 4.3 mm	37479 RP



Handle for Impl Rescue Collar & Guides 37480



### Trephine Drills

	Conical Connection	Tri-channel	External Hex
Ø 3.2/4.0 mm	37477 3.0		
Ø 3.8/4.6 mm	37928 NP	NP	NP
Ø 4.4/5.2 mm	37929 RP	RP, NobelReplace Platform Shift NP	RP
Ø 5.2/6.2 mm	37930 RP	WP, NobelReplace Platform Shift RP	WP
Ø 5.6/6.6 mm	37931 WP		
Ø 6.2/7.0 mm	37932	6.0, NobelReplace Platform Shift WP	WP

STERILE R



### Kits

Implant Retrieval Kit	37470
Implant Retrieval Kit Box	37514



## Abutment screw retrieval

### Abutment Screw Remover

	Conical Connection	Tri-channel	External Hex
3.0	37505		
NP	37506	37506	37506
RP, WP, 6.0	37507	37507	37507

STERILE R



### Rescue Drill Guide

	Conical Connection	Tri-channel	External Hex
3.0	37485		
NP	37486	37488	37482
RP	37487	37489	37483
WP	37933	37490	37484
6.0		37491	



### Abutment Screw Retrieval Instrument

	Conical Connection	Tri-channel	External Hex
3.0, NP	37503	37503	37503
RP, WP, 6.0	37504	37504	37504

STERILE R



### Screw Tap Repair

	Conical Connection	Tri-channel	External Hex
M1.4	37497 3.0		
M1.6	37498 NP		NP
M1.8	37499	NP	
M2	37517 RP, WP	RP, WP, 6.0	RP
M2.5	37500		WP

STERILE R



### Abutment Screw Retrieval Reverse Drill

	Conical Connection	Tri-channel	External Hex
3.0, NP	37501	37501	37501
RP, WP, 6.0	37502	37502	37502

STERILE R




Handle for Machine Instruments\* 29161



### Kits

Abutment Screw Retrieval Kit	37481
Abutment Screw Retrieval Kit Box	37515



## Abutment retrieval

### Abutment Retrieval Instrument

	Conical Connection
NP	Zr 37512
RP, WP	Zr 37882
NP	Ti 36247
RP, WP	Ti 37881



### Abutment Release Pin

	Conical Connection
3.0	37509
NP	37510
RP, WP	37511



### Kits

Abutment Retrieval Kit	37508
Abutment Retrieval Kit Box	37516



nobelbiocare.com/rescue