



**Goldstein
Crown Remover
Straight**
GCR0

For anterior
crown removal.



**Goldstein
Crown Remover
45° Angle**
GCR45

For cuspids,
bicusps,
and even
first molars.




**Goldstein
Crown Remover
Occlusal**
GCROS

For occlusal
separation
especially in
hard-to-remove
crowns that have
been bonded
to the tooth.



For permanent removal of crowns by breaking the seal between tooth and crown after sectioning with a bur. The special right angle handles are designed to torque the crown itself instead of destructive forces typically applied to the tooth which can lead to fracture.



**Goldstein
Crown Remover
Right Angle**
GCR90

For molars.

